

AEON™

RELOAD SELECTION

WIDEST
RANGE
OF
ANVILS



REGULAR TIP



SHORT TIP



CURVED TIP

RELOAD STAPLE HEIGHT

GRAY		WHITE		ORANGE		PURPLE		BLACK	
OPEN	CLOSED	OPEN	CLOSED	OPEN	CLOSED	OPEN	CLOSED	OPEN	CLOSED
2.0	0.75	2.5	1.0	3.25	1.5	4.0	1.8	5.0	2.2





30mm RELOADS



45mm RELOADS



60mm RELOADS

Reload Type	Cartridge Lengths	Product Codes*	Open Height	Closed Height
Gray	30, 45mm	AESR30G, AESR45G	2.0mm	0.75mm
Gray Curved Tip	30, 45mm	ASR30GC, AESC45G	2.0mm	0.75mm
Gray Short Tip	30mm	ASR30GS	2.0mm	0.75mm
White	30, 45, 60mm	AESR30W, AESR45W, AESR60W	2.5mm	1.0mm
White Curved Tip	30, 45, 60mm	ASR30WC, AESC45W, AESC60W	2.5mm	1.0mm
White Short Tip	30, 60mm	ASR30WS, ASR60WS	2.5mm	1.0mm
Orange	30, 45, 60mm	AESR30R, AESR45R, AESR60R	3.25mm	1.5mm
Orange Curved Tip	30, 45, 60mm	ASR30RC, AESC45R, AESC60R	3.25mm	1.5mm
Orange Short Tip	30, 60mm	ASR30RS, ASR60RS	3.25mm	1.5mm
Purple	30, 45, 60mm	AESR30P, AESR45P, AESR60P	4.0mm	1.8mm
Purple Curved Tip	30, 45, 60mm	ASR30PC, AESC45P, AESC60P	4.0mm	1.8mm
Purple Short Tip	30, 60mm	ASR30PS, ASR60PS	4.0mm	1.8mm
Black	45, 60mm	AESR45B, AESR60B	5.0mm	2.2mm
Black Curved Tip	45mm	ASR45BC	5.0mm	2.2mm

Handle Type	Shaft Length	Product Codes*
Short	60mm	AESH060, AETH060
Medium	160mm	AESH160, AETH160
Long	260mm	AESH260, AETH260

