

**FUJIFILM**

*FUJIFILM Electronic Video Bronchoscope*

**EB-530P**




Fujifilm provides for a new level of small caliber bronchoscopy. An impressive 3.8mm distal end diameter, combined with the 1.2mm forceps channel, allows for improved insertion capabilities into the peripheral bronchi.

Additionally, as with all Fujifilm bronchoscopes, superior imaging capabilities and a wide 120 degree field of view are combined to provide improved visual capabilities.



### Specifications

Field of view	120°
Observation range	3~100mm
Distal end diameter	φ 3.8mm
Flexible portion diameter	φ 3.8mm
Bending capability	UP/DOWN 180° /130°
Working length	600mm
Total length	890mm
Forceps channel diameter	φ 1.2mm



Generic Name: Bronchoscope, flexible, video